**EN INSTRUCTIONS FOR USE**

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# INSTRUCTIONS FOR USE

## Pedal exerciser

Product code: TIMgym


### 1.Characteristics

The lower limb pedal exerciser is made of a powder-coated steel frame. It is equipped with resistance mechanism, four anti-slip pads and foot straps. It has a rotary knob with a wide range of resistance adjustment, making it suitable for various patient groups. It can be used in hospitals, hospices and at home as a rehabilitation and personal recreation device.

#### Elements of the pedal exerciser:

- a. Pedal exerciser legs (2 pcs.)
- b. Pedal exerciser frame (1 pc.)
- c. Excerciser main mechanism (1 pc.)
- d. Adjustment knob (1 pc.)
- e. Screw with nuts type A (length 20 mm, 2 pcs.)
- f. Screws type B (length 10 mm, 4 pcs.)
- g. Screws washers (6 pcs.)

### 2. Application

The pedal exerciser is designed for patients with contracture or muscular atrophy caused by long-term immobilisation due to a plaster bandage or surgery, as well as for patients who need to restore motor functions and re-learn the correct movement pattern due to damage to the neurological system. It is possible to adjust the exerciser to the state of limb fitness and to the progress in rehabilitation owing to resistance adjustment.

### 3. Assembly and adjustment

- Insert the pedal exerciser leg into the frame and screw it with 2x B-type screws (placing washers beforehand). Screw the other pedal exerciser leg in a similar way.

**ATTENTION:** do not tighten the screws all the way.

- Once the legs are mounted to the frame, lightly press the frame to the ground and tighten all 4 mounting screws firmly. This will allow the pedal exerciser legs to level properly.
- Place the main mechanism on the pedal exerciser frame and fasten it with 2x screws and A-type nuts (placing washers beforehand).
- Screw the adjustment knob into the main mechanism.

### Resistance adjustment

In order to increase fitness levels and improve muscle work and muscle strength gains, the training intensity and/or duration of exercise should be gradually increased.

- To increase the resistance, turn the adjustment knob clockwise.
- To reduce the resistance, turn the adjustment knob anticlockwise.
- When performing low impact exercises, the adjustment knob should be completely unscrewed, whereas when performing resistance exercises, the adjustment knob should be tightened to the desired load.

### Usage

- Place the pedal exerciser on the floor and make sure the anti-slip pads are firmly in place to prevent the pedal exerciser from moving.
- Exercises should be carried out in a seated position, the pedal exerciser should be positioned exactly in front of the user so that the hip joint is bent at an angle of no less than 90°.

- Place feet on the pedals and adjust the strap to foot circumference, this will prevent the foot from slipping off the pedal during exercise.
- Set the resistance using the adjustment knob.

#### 4. Technical data

Base dimensions	40 cm x 50 cm
Base height from the ground	27 cm
Weight	2,35 kg
Maximum capacity	120 kg
Frame colour	black

#### 5. Notes

**Timago International Group does not bear responsibility of improper use of the pedal exerciser, failure to comply safety regulations and misuse.**

#### 6. Maintenance

Carry out regular maintenance of the technical condition of the pedal exerciser:

- It is recommended to clean the pedal exerciser regularly with a soft cloth and soapy water. After cleaning, dry the pedal exerciser thoroughly with soft cloth.
- Do not use strong polishing agents for cleaning as it may cause damage to the protective coatings.
- Remove impurities from moving elements of the exerciser on an ongoing basis.

#### 7. Storage and transport conditions

Device should not be exposed to direct sunlight, water or moisture.

## 8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

## 9. Labels

	Reference number		Manufacturer
	LOT number		Manufacture date
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

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