

## Walk Assist Jacket

- Designed to assist an ambulant patient with an unstable gait
- The jacket has soft handles positioned to enable the carers to offer some postural support and avoid grabbing clothing/limbs
- Carers can assist from behind or at the sides
- Ideal to assist those using walking sticks/frame
- Easy to apply like a jacket with zip fastening and side zips to allow expansion for outer clothing
- Mesh wash bag supplied for use when laundering jacket
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



## SIL-WALK



### Product Codes & Specifications:



|     | 100% Polyester Mesh | A           | B         | C    | D    |
|-----|---------------------|-------------|-----------|------|------|
| XXS | SIL-WALK-XXSM       | 59 - 77cm   | 38 - 44cm | 41cm | 53cm |
| XS  | SIL-WALK-XSM        | 77 - 87cm   | 44 - 50cm | 45cm | 58cm |
| S   | SIL-WALK-SM         | 87 - 107cm  | 50 - 57cm | 49cm | 62cm |
| M   | SIL-WALK-MED        | 97 - 117cm  | 57 - 64cm | 53cm | 72cm |
| L   | SIL-WALK-LRG        | 112 - 132cm | 64 - 69cm | 59cm | 83cm |

